





### Analysis of Bevacizumab and Erbitux Monoclonal Antibodies on Zenix™ SEC-300 (7.8 x 300 mm)



Column:	Zenix™ SEC-300 (3 μm, 300 Å, 7.8x300 mm)
Mobile Phase:	150 mM phosphate buffer, pH 7.0
Flow Rate:	1.0 mL/min;
Temperature:	25 °C
Detection:	UV 280 nm
Injection Volume:	10 μL
Samples:	Bevacizumab (trade name Avastin), a humanized monoclonal antibody, 150 kD



Column:	Zenix™ SEC-300 (3 μm, 300 Å, 7.8x300 mm)
Mobile Phase:	20 mM sodium phosphate buffer(pH7.0), 0.3 M NaCl
Flow Rate:	1.0 mL/min;
Temperature:	25 °C
Detection:	UV 214 nm
Injection Volume:	10 μL
Samples:	Erbitux (humanized monoclonal antibody)

Keywords: Size exclusion, Zenix, monoclonal antibody, bevacizumab, erbitux